

# Bacterial Transformation (activity)

## Contents

- [1 Bacterial Transformation](#)
- [2 Exercise: Transformation of Bacteria with RE Identified Plasmids](#)
- [3 Hypothesize: What will I expect of my transformed cells?](#)

## Bacterial Transformation

*Escherichia coli* are commensal gram negative bacteria found in the guts of humans. They have the capacity to double every twenty minutes and make a favorable carrier of recombinant DNA. Plasmid DNA can be introduced into *E. coli* easily after making them competent. One method to achieve this is through chemical competence with heat shock. In this process, the bacteria are incubated in  $\text{CaCl}_2$  solution on ice. The cold serves to slow down molecular motion of the plasma membrane while the  $\text{Ca}^{2+}$  ions remove the charge-charge repulsion between the phospholipids and the negatively charged DNA seeking to gain entry into the cell. Cells are placed for a short period of time at  $42^\circ\text{C}$  to induce **heat shock**. This heat shock results in the cell taking up the DNA. This method is very low efficiency so many bacteria do not take in any DNA. Cells are allowed to recover from heat shock at  $37^\circ\text{C}$  in rich nutrient broth to allow for the production of the antibiotic resistance proteins encoded on the vector as a selection marker. Transformed cells are then spread across an agar plate containing the antibiotic which will then kill all non-transformed cells. Only the bacteria containing the vector with the antibiotic resistance gene will survive and replicate to form small colonies on the surface of the agar.

- <http://www.dnalc.org/view/15918-Transformation.html>
- <http://www.dnalc.org/view/15916-DNA-transformation.html>

## Exercise: Transformation of Bacteria with RE Identified Plasmids

1. Each group retrieves the 2 miniprep plasmids from the previous week in the freezer and allow to thaw on ice.
2. Bring 2 agar plates to room temperature


2 plates will contain antibiotic, X-Gal and arabinose

3. For each plasmid, obtain 250  $\mu\text{l}$  of transformation buffer (50mM  $\text{CaCl}_2$ ) in microfuge tubes and place on ice for 10 minutes
4. Take an inoculating loop and remove a single colony of bacteria from a freshly streaked plate grown overnight
5. Swirl bacteria in each tube containing transforming solution to distribute bacteria throughout solution
6. Pipette 5  $\mu\text{l}$  of plasmid into the tube and incubate on ice for 10 minutes
7. During this incubation, flip the warmed plates and label them with your group names.


8. Place transformation tubes into 42°C heatblock for 1 minute to heat shock the cells
9. Add 500µl fresh SOC media (or LB) and incubate at 37°C for 15 minutes.
10. Pipette 150µl of transformation solution onto each plate and spread across the plate.
11. Turn plates agar side up and place them into 37°C incubator overnight. (your instructor will retrieve them and place them into refrigerator)

## Hypothesize: What will I expect of my transformed cells?


From the previous lab, we can identify our plasmids. The plasmids are either pGlo, pUC18/19 or pUC18/19 with a 6kb insert disrupting the LacZ gene. pGlo contains a gene that encodes the protein GFP that will fluoresce green under UV light and is 5.4kb. pUC is typically 2.7kb in size. LacZ is a gene encoding the protein  $\beta$ -Galactosidase, the enzyme that hydrolyzes lactose into the monosaccharides galactose and glucose. X-Gal is a chemical resembling lactose, however upon hydrolysis, the molecule deposits a blue coloring into the cell.



- [pUC19 genbank file](#)




- [pGlo genbank file](#)



- [pUC with insert file](#)

1. If the previously mentioned plasmids were digested by *EcoRI*, label the lanes below with the appropriate plasmid ([pGlo](#), [pUC](#), [pUC-inserted](#))
2. Predict if your transformants will be green under UV, white in all conditions or blue.
3. For additional help on this problem, utilize the [In silico digestion activity](#)



Tags: [visual communication](#)