

# 3 EXERCISE

## REVIEW SHEET


# The Microscope

Instructors may assign a portion of the Review Sheet questions using **Mastering A&P™**

Name FREEDOM ALMANZAR Lab Time/Date 3/6/21

### Care and Structure of the Compound Microscope

1. Label all indicated parts of the microscope.



2. Explain the proper technique for transporting the microscope.

YOU MUST HOLD IT TIGHT AND UPRIGHT POSITION, ONE HAND ON THE ARM AND OTHER ON BASE.

3. Each of the following statements is either true or false. If true, write *T* on the answer blank. If false, correct the statement by writing on the blank the proper word or phrase to replace the one that is underlined.

- F- SPECIAL GRIT-FREE PAPER 1. The microscope lens may be cleaned with any soft tissue.
- F- LOWEST POWER 2. The microscope should be stored with the oil immersion lens in position over the stage.
- T 3. When beginning to focus, use the scanning objective lens.
- AWAY FROM 4. When focusing on high power, always use the coarse adjustment knob to focus.
- T 5. A coverslip should always be used with wet mounts.

4. Match the microscope structures in column B with the statements in column A that identify or describe them.

**Column A**

- i 1. platform on which the slide rests for viewing
- D 2. used to adjust the amount of light passing through the specimen
- E 3. controls the movement of the slide on the stage
- B 4. delivers a concentrated beam of light to the specimen
- C 5. used for precise focusing once initial focusing has been done
- F 6. carries the objective lenses; rotates so that the different objective lenses can be brought into position over the specimen.

**Column B**

- a. coarse adjustment knob
- b. condenser
- c. fine adjustment knob
- d. iris diaphragm lever
- e. mechanical stage
- f. nosepiece
- g. objective lenses
- h. ocular lens
- i. stage

5. Define the following terms.

total magnification: BEING ABLE TO VIEW EQUALLY TO THE POWER OF THE OCULAR LENS MULTIPLIED BY THE POWER OF OBJECTIVE LENS USED.

resolution: ABILITY TO DISCRIMINATE TWO CLOSE OBJECT AS SEPARATE.

### Viewing Objects Through the Microscope

6. Complete, or respond to, the following statements:

- WORKING DISTANCE 1. The distance from the bottom of the objective lens to the surface of the slide is called the \_\_\_\_\_.
- IN THE LEFT 2. Assume there is an object on the left side of the field that you want to bring to the center (that is, toward the apparent right). In what direction would you move your slide? \_\_\_\_\_.
- FIELD 3. The area of the slide seen when looking through the microscope is the \_\_\_\_\_.
- 93 4. If a microscope has a 10× ocular lens and the total magnification is 950×, the objective lens in use at that time is \_\_\_\_\_ ×.


INCREASE CONTRAST 5. Why should the light be dimmed when looking at living (nearly transparent) cells?

PARFOCAL 6. If, after focusing in low power, you need to use only the fine adjustment to focus the specimen at the higher powers, the microscope is said to be \_\_\_\_\_.

3 7. You are using a 10× ocular and a 15× objective, and the field diameter is 1.5 mm. The approximate field size with a 30× objective is \_\_\_\_\_ mm.

.4 8. If the diameter of the low-power field is 1.5 mm, an object that occupies approximately a third of that field has an estimated diameter of \_\_\_\_\_ mm.

7. You have been asked to prepare a slide with the letter *F* on it (as shown below). In the circle below, draw the *F* as seen in the low-power field.




8. Estimate the length (longest dimension) of the object in μm:

$$(1.6 \text{ mm} \times 100) / 100 =$$

Total magnification = 100×

Field diameter = 1.6 mm

Length of object = 1.6 μm



9. Say you are observing an object in the low-power field. When you switch to high power, it is no longer in your field of view.

Why might this occur? FIELD DECREASES AS MAGNIFICATION INCREASES. THE OBJECT TO BE OUTSIDE IS ALSO HAS HIGHER FIELD

What should you do initially to prevent this from happening? CENTER THE OBJECT

10. Do the following factors increase or decrease as one moves to higher magnifications with the microscope?

resolution: INCREASE amount of light needed: INCREASE

working distance: DECREASE depth of field: DECREASE

11. A student has the high-power lens in position and appears to be intently observing the specimen. The instructor, noting a working distance of about 1 cm, knows the student isn't actually seeing the specimen.

How so? WORKING DISTANCE FOR HIGH POWER LENS IS 1MM

12. Describe the proper procedure for preparing a wet mount.

PLACE SPECIMEN ON THE SLIDE AND PLACE A DROP OF WATER OR SALINE IN THE SLIDE, MIX THE SPECIMEN USING A TOOTHPICK. THEN ADD A DROP OF STAIN AND MIX CAREFULLY LOWER THE COVERSLIP ONTO THE SPECIMEN

13. Indicate the probable cause of the following situations during use of a microscope.

a. Only half of the field is illuminated: THE LENS IS NOT CORRECTLY ROTATED INTO PLACE.

b. The visible field does not change as the mechanical stage is moved: THE SLIDE IS NOT CORRECTLY POSITIONED IN THE CLAMP

14. **+** A blood smear is used to diagnose malaria. In patients with malaria, the protozoa can be found near and inside red blood cells. Explain why a microscope capable of high magnification and high resolution would be needed to diagnose malaria.

TO BE ABLE TO SEE INSIDE THE RED BLOOD CELL

15. **+** Histopathology is the use of microscopes to view tissues to diagnose and track the progression of diseases. Why are thin slices of tissue ideal for this procedure?

IT'S IDEAL BECAUSE THE LIGHT IS ABLE TO PASS THROUGH THE SAMPLE