

Brianna Gorham

Professor Haque

Lab

Investigating Diffusion and Osmosis Through Nonliving Membranes

Materials:

- 2 dialysis sacs containing 20 ml of 40% glucose solution
- 1 dialysis sac containing 20 ml of 10% NaCl solution
- 1 dialysis sac containing 20 ml of 40% sucrose solution
- 3 beakers $\frac{1}{2}$ filled with distilled water
- 1 beaker $\frac{1}{2}$ filled with 40% glucose solution
- Boiling water
- Benedict's solution (light blue color)
- Test tubes
- Small funnel
- Fine Twine
- 25-ml graduated cylinder
- Wax marking pencil
- Laboratory Balance
- Test tube rack

Procedure/ Instructions

- Label all materials with wax marker pencil. (All sacs label as 1A-4A and all beakers as 1B-4B)
- Using a funnel, Fill each sac one at a time, as the following: Sac 1A= fill with 20 ml of 20% of glucose solution Sac 2B= fill with 20ml of 40% glucose solution Sac 3C= fill with 20ml of 10% of NaCl solution Sac 4D= fill with 20 ml of 40% of sucrose solution

- Tie each sack with fine twine
- Weigh each sac with a laboratory balance to make sure everything is perfectly equipped; Record weight accordingly in the chart.
- Fill each beaker $\frac{1}{2}$ as the following: Beaker 1B= $\frac{1}{2}$ of distilled water Beaker 2B= $\frac{1}{2}$ of 40% of glucose solution Beaker 3B= $\frac{1}{2}$ of distilled water Beaker 4B= $\frac{1}{2}$ of distilled water
- Drop the sac into the appropriate beaker, making sure sac is completely covered by the beaker solution.
- Allow all sacs to remain still in beakers for 45 minutes; Observe.
- One by one, quickly, and carefully remove sac from beaker, blot sac dry and record final weight accordingly in the data chart.
- Calculate the weight change (increase/decrease, by how much) for each sac.
- After 45 minutes, boil a beaker of water; Obtain dropper bottles of Benedicts solution, silver nitrate solution, test tube rack, 4 test tubes, and a test tube holder.
- Now we will test for sugar: Obtain sample (will show as colorless) and Benedicts Solution (will appear as light blue color) and mix. Only test #1,2, and 4
- Place the test tube in the boiling water
- Observe the color change (will appear yellow, green, or brown)
- Record down your observations
- Test #3 testing for NaCl
- Record all observations and conclusions

Results:

Beaker	Contents Of Sac	Initial Weight	Final Weight	Weight Change	Test-Beaker Fluid	Test-Fluid Sac Fluid	
Beaker 1 <i>½ filled with distilled water</i>	Sac 1, 20 ml of 40% glucose solution	7.1 gm	8.0gm	Increased 0.9g	Benedict's Test: Positive	Benedict's Test : Positive	Positive
Beaker 2 <i>½ filled with 40% glucose solution</i>	Sac 2, 20 ml of 40% glucose solution	6.9gm	6.9gm	No Change	Positive	Positive	Positive
Beaker 3 <i>½ filled with distilled water</i>	Sac 3, 20 ml of 10% NACI Solution	7.2gm	7.8gm	Increased 0.6g	AgNO3 test N/A	N/A	Positive
Beaker 4 <i>½ filled with distilled water</i>	Sac 4, 20 ml of 40% sucrose solution	7.1gm	8.0gm	Increased 0.9	Benedict's Test:	Negative	Positive

Hypothesis: If the substance is added to the dialysis bag and is placed into the water then the water will seep into the dialysis bag will have a great mass this process is called diffusion